

For Research Use Only

NEO1/Neogenin Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88292-2-PBS

Basic Information

Catalog Number: 88292-2-PBS	GenBank Accession Number: NM_002499.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 4756	CloneNo.: 260336E2
Isotype: IgG	UNIPROT ID: Q92859-1	
Immunogen Catalog Number: EG5518	Full Name: neogenin homolog 1 (chicken)	
	Calculated MW: 160 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Neogenin concentration was determined to be 22.19 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 20.1 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load, 0.2 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03454-1, Neogenin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88292-3-PBS. Detection antibody: 88292-2-PBS. Standard: Eg5518. Range: 0.312-40 ng/mL.