

For Research Use Only

CPA1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88295-3-PBS



Basic Information

Catalog Number:

88295-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8586

GenBank Accession Number:

BC005279

GeneID (NCBI):

1357

UNIPROT ID:

P15085

Full Name:

carboxypeptidase A1 (pancreatic)

Calculated MW:

419 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

260663D2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of six individual healthy human donors was measured. The CPA1 concentration of detected samples was determined to be 5.51 ng/mL with a range of 0.70-17.97 ng/mL.



Sandwich ELISA standard curve of MP03531-1, Human CPA1 Recombinant Matched Antibody Pair - PBS only. 88295-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag8586. 88295-2-PBS was HRP conjugated as the detection antibody. Range: 0.977-62.5 ng/mL.