

For Research Use Only

DMKN Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88298-2-PBS



Basic Information

Catalog Number:

88298-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9242

GenBank Accession Number:

BC011886

GeneID (NCBI):

93099

UNIPROT ID:

Q6E0U4

Full Name:

dermokine

Calculated MW:

476 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

260498B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DMKN concentration was determined to be 35.2 ng/mL in SK-MEL-30 cell extract based on a 1.2 mg/mL extract load, 2.1 ng/mL in MCF-10A cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03499-1, DMKN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88298-3-PBS. Detection antibody: 88298-2-PBS. Standard: Ag9242. Range: 0.781-100 ng/mL.