

For Research Use Only

Sam68 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88301-3-PBS



Basic Information

Catalog Number:

88301-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0235

GenBank Accession Number:

BC000717

GeneID (NCBI):

10657

UNIPROT ID:

Q07666

Full Name:

KH domain containing, RNA binding, signal transduction associated 1

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

260478D5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Sam68 concentration was determined to be 17.46 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, 6.52 ng/mL in Caco-2 cell extract based on a 1.2 mg/mL extract load and 683.88 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03453-2, Human Sam68 Recombinant Matched Antibody Pair - PBS only. 88301-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0235. 88301-3-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL