

For Research Use Only

HSF1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88318-1-PBS



Basic Information

Catalog Number: 88318-1-PBS	GenBank Accession Number: BC014638	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3297	CloneNo.: 260631H4
Isotype: IgG	UNIPROT ID: Q00613	
Immunogen Catalog Number: AG0739	Full Name: heat shock transcription factor 1	
	Calculated MW: 529 aa, 57 kDa	
	Observed MW: 68-80 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

HSF1 belongs to heat-shock transcription factors that activate heat-shock response genes under conditions of heat or other stresses. Also, HSF1 has been linked with oogenesis, spermatogenesis, and placental development. It can activate AKT and inactivate JNK and CASP3 to protect cardiomyocytes from death. And it has a role in the regulation of life span and establishes a role for SIRT1 in protein homeostasis and heat-shock response. The calculated molecular weight of HSF1 is 57 kDa, but HSF1 migrates at approximately 80 kDa which likely represents different phosphorylation states (PMID: 18434628).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88318-1-RR (HSF1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88318-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded Human lung cancer (Adenocarcinoma) slide using 88318-1-RR (HSF1 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 88318-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HSF1 antibody (88318-1-RR, Clone: 260631H4) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 88318-1-PBS in a different storage buffer formulation.