

For Research Use Only

PRPF19 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88331-3-PBS



Basic Information

Catalog Number:	88331-3-PBS	GenBank Accession Number:	BC008719	Purification Method:	Protein A purification
Source:	Rabbit	GeneID (NCBI):	27339	CloneNo.:	260837B4
Isotype:	IgG	UNIPROT ID:	Q9UMS4		
Immunogen Catalog Number:	AG7700	Full Name:	PRP19/PSO4 pre-mRNA processing factor 19 homolog (S. cerevisiae)		
		Calculated MW:	55 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PRPF19 concentration was determined to be 229.4 ng/mL in HL-60 cell extract based on a 1.2 mg/mL extract load, 75.8 ng/mL in SKOV-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03457-1, PRPF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88331-3-PBS. Detection antibody: 88331-2-PBS. Standard: Ag7700. Range: 1.875-60 ng/mL.