

For Research Use Only

LSR Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88350-3-PBS



Basic Information

Catalog Number:

88350-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12961

GenBank Accession Number:

BC004381

GeneID (NCBI):

51599

UNIPROT ID:

Q86X29

Full Name:

lipolysis stimulated lipoprotein
receptor

Calculated MW:

71 kDa

Purification Method:

Protein A purification

CloneNo.:

260839B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LSR concentration was determined to be 1600.1 ng/mL in HT-29 cell extract based on a 1.2 mg/mL extract load, 852.5 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, THP-1 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03491-1, LSR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88350-3-PBS. Detection antibody: 88350-2-PBS. Standard: Ag12961. Range: 0.781-100 ng/mL.