

For Research Use Only

LSR Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88350-4-PBS



Basic Information

Catalog Number:

88350-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12961

GenBank Accession Number:

BC004381

GeneID (NCBI):

51599

UNIPROT ID:

Q86X29

Full Name:

lipolysis stimulated lipoprotein
receptor

Calculated MW:

71 kDa

Purification Method:

Protein A purification

CloneNo.:

260839D4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LSR concentration was determined to be 116.18 ng/mL in A431 cell extract based on a 1.20 mg/mL extract load, 6.22 ng/mL in THP-1 cell extract based on a 1.20 mg/mL extract load and 4.82 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP03491-2, Human LSR Recombinant Matched Antibody Pair - PBS only. 88350-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag12961. 88350-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.