

For Research Use Only

ITIH4 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88354-4-PBS



Basic Information

Catalog Number:

88354-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG10205

GenBank Accession Number:

NM_002218.4

GeneID (NCBI):

3700

UNIPROT ID:

Q14624-1

Full Name:

inter-alpha (globulin) inhibitor H4
(plasma Kallikrein-sensitive
glycoprotein)

Calculated MW:

103 kDa

Purification Method:

Protein A purification

CloneNo.:

260966D1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of six individual healthy human donors was measured. The ITIH4 concentration of detected samples was determined to be 43.0 ug/mL with a range of 27.7-62.4 ug/mL.



Sandwich ELISA standard curve of MP03496-2, Human ITIH4 Recombinant Matched Antibody Pair - PBS only. 88354-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg10205. 88354-3-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.