

For Research Use Only

GOLSYN Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88360-2-PBS



Basic Information

Catalog Number:

88360-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9344

GenBank Accession Number:

BC012082

GeneID (NCBI):

55638

UNIPROT ID:

Q9NX95

Full Name:

Golgi-localized protein

Calculated MW:

662 aa, 72 kDa

Purification Method:

Protein A purification

CloneNo.:

260730A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Syntabulin concentration was determined to be 4.8 ng/mL in HT-29 cell extract based on a 1.2 mg/mL extract load, 2.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, HL-60 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03456-1, Syntabulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88360-2-PBS. Detection antibody: 88360-1-PBS. Standard: Ag9344. Range: 0.312-40 ng/mL.