

For Research Use Only

EIF4B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88367-2-PBS



Basic Information

Catalog Number:

88367-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12325

GenBank Accession Number:

BC073154

GeneID (NCBI):

1975

UNIPROT ID:

P23588

Full Name:

eukaryotic translation initiation factor 4B

Calculated MW:

611 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

260892B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean EIF4B concentration was determined to be 247.9 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load, 122.5 ng/mL in HEK-293 cell extract based on a 1.0 mg/mL extract load, 87.1 ng/mL in MCF-7 cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03486-1, EIF4B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88367-3-PBS. Detection antibody: 88367-2-PBS. Standard: Ag12325. Range: 0.156-20 ng/mL.