

For Research Use Only

# ZCCHC6 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88391-1-PBS



## Basic Information

Catalog Number:

88391-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18240

GenBank Accession Number:

BC032456

GeneID (NCBI):

79670

UNIPROT ID:

Q5VYS8

Full Name:

zinc finger, CCHC domain containing 6

Calculated MW:

1495 aa, 171 kDa

Observed MW:

171 kDa

Purification Method:

Protein A purification

CloneNo.:

260850G10

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

ZCCHC6 also named as KIAA1711, HS2 or TUT7 is a 1495 amino acid protein, which contains three CCHC-type zinc fingers and two PAP-associated domains. A member of the DNA polymerase type-B-like family, ZCCHC6 as a uridylyltransferase that mediates RNA uridylation.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88391-1-RR (ZCCHC6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88391-1-PBS in a different storage buffer formulation.