

For Research Use Only

FAF2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88402-2-PBS



Basic Information

Catalog Number:

88402-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33251

GenBank Accession Number:

BC014001

GeneID (NCBI):

23197

UNIPROT ID:

Q96CS3

Full Name:

Fas associated factor family member

2

Calculated MW:

445 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

260728F10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FAF2 concentration was determined to be 35.2 ng/mL in HepG2 cell extract based on a 1.0 mg/mL extract load, 59.2 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load, 8.1 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03535-1, FAF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88402-3-PBS. Detection antibody: 88402-2-PBS. Standard: Ag33251. Range: 0.781-100 ng/mL.