

For Research Use Only

AKAP7 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88409-2-PBS



Basic Information

Catalog Number:

88409-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3297

GenBank Accession Number:

BC016927

GeneID (NCBI):

9465

UNIPROT ID:

O43687

Full Name:

A kinase (PRKA) anchor protein 7

Calculated MW:

81 aa, 9 kDa

Observed MW:

9-15 kDa

Purification Method:

Protein A purification

CloneNo.:

260933E6

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

AKAP7, also named as AKAP15 and AKAP18, localize PKA in complexes with CaV1.2 and PLN, respectively. AKAP7 targets the cAMP-dependent protein kinase (PKA) to the plasma membrane, and permits functional coupling to the L-type calcium channel. AKAP7 has short isoform AKAP7 alpha and AKAP7 beta with MW 9-15 kDa; long isoform AKAP7 gamma with MW 37-42 kDa. AKAP7alpha was highly expressed only in brain and weakly in lung lysates from WT animals.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88409-2-RR (AKAP7 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88409-2-PBS in a different storage buffer formulation.



IP result of anti-AKAP7 (IP:88409-2-RR, 4ug; Detection:88409-2-RR 1:1000) with human brain tissue lysate 1200 ug. This data was developed using the same antibody clone with 88409-2-PBS in a different storage buffer formulation.