

For Research Use Only

AKAP7 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88409-3-PBS



Basic Information

Catalog Number: 88409-3-PBS	GenBank Accession Number: BC016927	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 9465	CloneNo.: 260933F1
Isotype: IgG	UNIPROT ID: O43687	
Immunogen Catalog Number: AG3297	Full Name: A kinase (PRKA) anchor protein 7	
	Calculated MW: 81 aa, 9 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean AKAP7 concentration was determined to be 60.7 ng/mL in RT-4 cell extract based on a 1.2 mg/mL extract load, HaCaT cell extract were served as a negative control.



Cytometric bead array standard curve of MP03506-1, AKAP7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88409-1-PBS. Detection antibody: 88409-3-PBS. Standard: Ag3297. Range: 0.156-20 ng/mL.