

For Research Use Only

UNC93B1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88417-3-PBS



Basic Information

Catalog Number:

88417-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27937

GenBank Accession Number:

BC101568

GeneID (NCBI):

81622

UNIPROT ID:

Q9H1C4

Full Name:

unc-93 homolog B1 (C. elegans)

Calculated MW:

597 aa, 67 kDa

Purification Method:

Protein A purification

CloneNo.:

260829A11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean UNC93B1 concentration was determined to be 72.87 ng/mL in HL-60 cell extract based on a 1.20 mg/mL extract load and 8.59 ng/mL in IMR-32 cell extract based on a 1.70 mg/mL extract load.



Sandwich ELISA standard curve of MP03543-2, Human UNC93B1 Recombinant Matched Antibody Pair - PBS only. 88417-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27937. 88417-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.