

For Research Use Only

BLVRA Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88427-3-PBS



Basic Information

Catalog Number:

88427-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1166

GenBank Accession Number:

BC008456

GeneID (NCBI):

644

UNIPROT ID:

P53004

Full Name:

biliverdin reductase A

Calculated MW:

33.2 kDa

Purification Method:

Protein A purification

CloneNo.:

260945C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BLVRA concentration was determined to be 48.7 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load, 100.6 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, Ramos cell extract were served as a negative control.



Cytometric bead array standard curve of MP03540-1, BLVRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88427-3-PBS. Detection antibody: 88427-2-PBS. Standard: Ag1166. Range: 0.469-60 ng/mL.