

For Research Use Only

BHLHE41 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88440-4-PBS



Basic Information

Catalog Number:

88440-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3381

GenBank Accession Number:

BC025968

GeneID (NCBI):

79365

UNIPROT ID:

Q9C0J9

Full Name:

basic helix-loop-helix family,
member e41

Calculated MW:

482 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

260849C11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BHLHE41 concentration was determined to be 1,884.8 pg/mL in RT-4 cell extract based on a 1.2 mg/mL extract load, 673.3 pg/mL in U2OS cell extract based on a 1.2 mg/mL extract load and 1,126.3 pg/mL in A431 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03609-1, Human BHLHE41 Recombinant Matched Antibody Pair - PBS only. 88440-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3381. 88440-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.