

For Research Use Only

RAP1B Recombinant monoclonal antibody, PBS Only

Catalog Number: 88443-1-PBS



Basic Information

Catalog Number: 88443-1-PBS	GenBank Accession Number: BC000176	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5908	CloneNo.: 261087A12
Isotype: IgG	UNIPROT ID: P61224	
Immunogen Catalog Number: AG1282	Full Name: RAP1B, member of RAS oncogene family	
	Calculated MW: 21 kDa	
	Observed MW: 21 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

RAP1B is a member of the RAS-like small GTP-binding protein superfamily. Members of this family regulate multiple cellular processes including cell adhesion and growth and differentiation. RAP1B localizes to cellular membranes and has been shown to regulate integrin-mediated cell signaling. It also plays a role in regulating outside-in signaling in platelets.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88443-1-RR (RAP1B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88443-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using RAP1B antibody (88443-1-RR, Clone: 261087A12) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 88443-1-PBS in a different storage buffer formulation.