

For Research Use Only

CENPA Recombinant monoclonal antibody, PBS Only

Catalog Number: 88445-1-PBS



Basic Information

Catalog Number:

88445-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25025

GenBank Accession Number:

BC002703

GeneID (NCBI):

1058

UNIPROT ID:

P49450

Full Name:

centromere protein A

Calculated MW:

16 kDa

Observed MW:

16-20 kDa

Purification Method:

Protein A purification

CloneNo.:

260865G5

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Centromere protein-A (CENP-A) is a Histone H3-like protein that contains a C-terminal H3-like domain, required for centromere localization of CENP-A, and an antigenic N-terminal domain. CENP-A, originally isolated from HeLa cells, is essential for kinetochore targeting of CENP-C. In the presence of DNA, CENP-A forms an octameric complex with histones H4, H2A and H2B. CENP-A specifically localizes to active centromeres and is a component of specialized centromeric nucleosomes, on which kinetochores are assembled. CENP-A is essential for nucleosomal packaging of centromeric DNA at interphase and functions as a centromere formation marker on the chromosome.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88445-1-RR (CENPA antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88445-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed PC-3 cells using CENPA antibody (88445-1-RR, Clone: 260865G5) at dilution of 1:200 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), Alpha Tubulin antibody (66031-1-Ig, Clone: 1E4C11, red). This data was developed using the same antibody clone with 88445-1-PBS in a different storage buffer formulation.