

For Research Use Only

UQCRFS1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88500-1-PBS



Basic Information

Catalog Number:

88500-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13315

GenBank Accession Number:

BC010035

GeneID (NCBI):

7386

UNIPROT ID:

P47985

Full Name:

ubiquinol-cytochrome c reductase,
Rieske iron-sulfur polypeptide 1

Calculated MW:

30 kDa

Observed MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

261156A1

Applications

Tested Applications:

WB, IF/ICC, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

The ubiquinol cytochrome c reductase UQCRFS1 is a key subunit of the cytochrome bc1 complex (complex III) of the mitochondrial respiratory chain. The UQCRFS1 gene encode a 30 kDa protein with a transit peptide and it appears to be involved in development of more aggressive phenotype of breast cancer (PMID:15047214). UQCRFS1 is also named as RISP, complex III subunit 5, cytochrome b-c1 complex subunit 5 and rieske iron-sulfur protein.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88500-1-RR (UQC RF S1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88500-1-PBS in a different storage buffer formulation.



IP result of anti-UQC RF S1 (IP:88500-1-RR, 4ug; Detection:88500-1-RR 1:3000) with HeLa cells lysate 1080 ug. This data was developed using the same antibody clone with 88500-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed MCF-7 cells using UQC RF S1 antibody (88500-1-RR, Clone: 261156A1) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 88500-1-PBS in a different storage buffer formulation.