

For Research Use Only

JUP Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88895-3-PBS



Basic Information

Catalog Number:

88895-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1627

GenBank Accession Number:

BC011865

GeneID (NCBI):

3728

UNIPROT ID:

P14923

Full Name:

junction plakoglobin

Calculated MW:

82 kDa

Purification Method:

Protein A purification

CloneNo.:

261780C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Gamma Catenin concentration was determined to be 145.6 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 76.2 ng/mL in T-47D cell extract based on a 1.2 mg/mL extract load, 1.5 ng/mL in A-204 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03731-1, Gamma Catenin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88895-3-PBS. Detection antibody: 88895-2-PBS. Standard: Ag1627. Range: 0.312-40 ng/mL.