

For Research Use Only

NFATC3 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88911-3-PBS



Basic Information

Catalog Number:

88911-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12985

GenBank Accession Number:

BC001050

GeneID (NCBI):

4775

UNIPROT ID:

Q12968

Full Name:

nuclear factor of activated T-cells,
cytoplasmic, calcineurin-dependent 3

Calculated MW:

116 kDa

Purification Method:

Protein A purification

CloneNo.:

261643H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NFATC3 concentration was determined to be 4.4 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 1.4 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03718-1, NFATC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88911-3-PBS. Detection antibody: 88911-2-PBS. Standard: Ag12985. Range: 0.078-10 ng/mL.