

For Research Use Only

FBXW7 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88954-3-PBS



Basic Information

Catalog Number:

88954-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28460

GenBank Accession Number:

NM_001349798

GeneID (NCBI):

55294

UNIPROT ID:

Q969H0

Full Name:

F-box and WD repeat domain containing 7

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

261793C9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FBXW7 concentration was determined to be 80.53 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load, 104.83 ng/mL in HEK-293 cell extract based on a 1.20 mg/mL extract load, 82.72 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load and 50.78 ng/mL in A431 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP03838-1, Human FBXW7 Recombinant Matched Antibody Pair - PBS only. 88954-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag28460. 88954-2-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.