

For Research Use Only

LMAN1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88996-2-PBS



Basic Information

Catalog Number:

88996-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4183

GenBank Accession Number:

BC032330

GeneID (NCBI):

3998

UNIPROT ID:

P49257

Full Name:

lectin, mannose-binding, 1

Calculated MW:

510 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

262111H6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ERGIC-53 concentration was determined to be 759.0 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 170.6 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03757-1, ERGIC-53 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88996-3-PBS. Detection antibody: 88996-2-PBS. Standard: Ag4183. Range: 0.312-40 ng/mL.