

For Research Use Only

KLF10 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 89007-2-PBS



Basic Information

Catalog Number:

89007-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31772

GenBank Accession Number:

BC011538

GeneID (NCBI):

7071

UNIPROT ID:

Q13118

Full Name:

Kruppel-like factor 10

Calculated MW:

469 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

261049C5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean KLF10 concentration was determined to be 1,998.95 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load and 655.48 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03712-1, Human KLF10 Recombinant Matched Antibody Pair - PBS only. 89007-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag31772. 89007-1-PBS was HRP conjugated as the detection antibody. Range: 3.13-100 ng/mL