

For Research Use Only

MRPS27 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 89044-4-PBS



Basic Information

Catalog Number:

89044-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11047

GenBank Accession Number:

BC030521

GeneID (NCBI):

23107

UNIPROT ID:

Q92552

Full Name:

mitochondrial ribosomal protein S27

Calculated MW:

168 aa, 19 kDa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

262016C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MRPS27 concentration was determined to be 7.2 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 2.7 ng/mL in THP-1 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03788-1, MRPS27 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 89044-4-PBS. Detection antibody: 89044-3-PBS. Standard: Ag11047. Range: 0.156-10 ng/mL.