

For Research Use Only

# Anti-Mouse CXCL1 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98047-1-PBS



## Basic Information

Catalog Number:

98047-1-PBS

Size:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008176

GeneID (NCBI):

14825

UNIPROT ID:

P12850-1

Full Name:

chemokine (C-X-C motif) ligand 1

Calculated MW:

10kd

Purification Method:

Protein A purification

CloneNo.:

240268D2

## Applications

Tested Applications:

FC (Intra)

Species Specificity:

mouse

## Background Information

CXCL1 is a member of CXC family and also known as keratinocyte-derived chemokines (KC) or growth-related oncogene (GRO). CXCL1 is expressed by macrophages, neutrophils and epithelial cells. CXCL1 binding specifically to the CXC chemokine receptor CXCR2, is involved in fibrogenesis and angiogenesis. This protein also plays a role in inflammation and as a chemoattractant for neutrophils. CXCL1 is upregulated in some types of human cancer, including colorectal, bladder, breast, prostate and skin cancers.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> untreated or LPS and Brefeldin A treated mouse splenocytes were intracellularly stained with 0.125 ug Anti-Mouse CXCL1 Rabbit Recombinant Antibody. Cells were then stained with APC Anti-Mouse TNF α . Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 98047-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98047-1-RR against Mouse CXCL1 were performed. The affinity constant is 35.4 nM.