

For Research Use Only

Anti-Mouse PD-1/CD279 Mouse Recombinant Antibody, PBS Only

Catalog Number: 98080-6-PBS



Basic Information

Catalog Number:

98080-6-PBS

Concentration:

1mg, 3 mg/ml

Source:

Mouse

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

18566

Full Name:

programmed cell death 1

Purification Method:

Protein A purification

CloneNo.:

29F.1A12

Applications

Tested Applications:

FC

Species Specificity:

mouse

Background Information

Programmed cell death 1 (PD-1, also known as CD279) is an immunoinhibitory receptor that belongs to the CD28/CTLA-4 subfamily of the Ig superfamily. It is a 288 amino acid (aa) type I transmembrane protein composed of one Ig superfamily domain, a stalk, a transmembrane domain, and an intracellular domain containing an immunoreceptor tyrosine-based inhibitory motif (ITIM) as well as an immunoreceptor tyrosine-based switch motif (ITSM) (PMID: 18173375). PD-1 is expressed during thymic development and is induced in a variety of hematopoietic cells in the periphery by antigen receptor signaling and cytokines (PMID: 20636820). Engagement of PD-1 by its ligands PD-L1 or PD-L2 transduces a signal that inhibits T-cell proliferation, cytokine production, and cytolytic function (PMID: 19426218). It is critical for the regulation of T cell function during immunity and tolerance. Blockade of PD-1 can overcome immune resistance and also has been shown to have antitumor activity (PMID: 22658127; 23169436). RMP1-14 is a monoclonal antibody that targets the murine PD-1 protein and can be used in mouse models to block the reaction between PD-1 and PD-L1/PD-L2 (PMID: 14764726; 31832706).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or anti-CD3/CD28 treated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse PD-1/CD279 Mouse RecAb (98080-6-MR, Clone: 29F.1A12) and PE-conjugated Goat Anti-Mouse IgG2a. Cells were then stained with CoraLite® Plus 647 Anti-Mouse CD3. Cells were not fixed. This data was developed using the same antibody clone with 98080-6-PBS in a different storage buffer formulation.

1x10⁶ anti-CD3/CD28 treated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse PD-1/CD279 Mouse RecAb (98080-6-MR, Clone: 29F.1A12) (red) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) (blue), and PE-conjugated Goat Anti-Mouse IgG2a. 1x10⁶ untreated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse PD-1/CD279 Mouse RecAb (98080-6-MR, Clone: 29F.1A12) (black, dashed) and PE-conjugated Goat Anti-Mouse IgG2a. Cells