

For Research Use Only

Anti-Human CD177 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98086-1-PBS



Basic Information

Catalog Number:

98086-1-PBS

Size:

1mg, 2mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC029167

GeneID (NCBI):

57126

UNIPROT ID:

Q8N6Q3

Full Name:

CD177 molecule

Calculated MW:

437 aa, 46 kDa

Purification Method:

Protein A purification

CloneNo.:

241339F4

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

CD177, also known as NB1, HNA-2a or PRV-1, is a glycosylphosphatidylinositol (GPI)-anchored glycoprotein that belongs to the uPAR/CD59/Ly6 snake toxin superfamily (PMID: 17898576). It is expressed exclusively by neutrophils, neutrophilic metamyelocytes, and myelocytes (PMID: 17580308). CD177 interacts with β 2 integrin CD11b/CD18 and CD11a/CD18, PECAM-1, and participates in neutrophil activation and transmigration (PMID: 17580308; 21193407). CD177 has been identified as a useful biomarker for myeloproliferative diseases (PMID: 17898576).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were surface stained with 0.25 ug Anti-Human CD177 Rabbit Recombinant Antibody (98086-1-RR, Clone:241339F4) or Isotype Control and PE-conjugated Goat Anti-Rabbit IgG. Cells were incubated with FC Receptor Block prior to fixing and staining. Cells were not fixed. This data was developed using the same antibody clone with 98086-1-PBS in a different storage buffer formulation.



1x10⁶ human peripheral blood leukocytes were surface stained with 0.25 ug Anti-Human CD177 Rabbit Recombinant Antibody (98086-1-RR, Clone:241339F4) (red) or Isotype Control (blue) and PE-conjugated Goat Anti-Rabbit IgG. Cells were incubated with FC Receptor Block prior to fixing and staining. Cells were not fixed. This data was developed using the same antibody clone with 98086-1-PBS in a different storage buffer formulation.