

For Research Use Only

Anti-Human CEACAM3/6 (CD66c/d) Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98169-1-PBS



Basic Information

Catalog Number:

98169-1-PBS

Size:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC106728

GeneID (NCBI):

1084

UNIPROT ID:

P40198

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 3

Calculated MW:

252 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

241254D1

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

Carcinoembryonic antigen-related cell adhesion molecules (CEACAMs) belong to the immunoglobulin superfamily of cell adhesion molecules. CEACAMs play major roles in cell-cell and cell-ECM adhesion and have been implicated in controlling cell proliferation, angiogenesis, and tissue remodeling (PMID: 23021083). CEACAM3, also named as CD66d and CGM1, is exclusively expressed on granulocytes and acts as the major granulocyte receptor mediating recognition and efficient opsonin-independent phagocytosis of CEACAM-binding microorganisms (PMID: 18606569). CEACAM6, also known as CD66c, is normally expressed on the surface of myeloid and epithelial surfaces and upregulated preferentially at the apical/luminal membranes of many tumors (PMID: 27595061; 35082925). This antibody is raised against 35-155 aa of human CEACAM3/CD66d and can cross-react with CEACAM6/CD66c.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were surface stained with 0.25 ug Anti-Human CEACAM3/6 (CD66c/d) Rabbit Recombinant Antibody (98169-1-RR, Clone: 241254D1) or 0.25 ug Isotype Control, and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed. This data was developed using the same antibody clone with 98169-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98169-1-RR against Human CEACAM3/6 (CD66c/d) were performed. The affinity constant is 2.19 nM.