

For Research Use Only

Anti-Mouse CD28 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98195-1-PBS



Basic Information

Catalog Number:

98195-1-PBS

Size:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007642.4

GeneID (NCBI):

12487

UNIPROT ID:

P31041

Full Name:

CD28 antigen

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

241707H6

Applications

Tested Applications:

FC

Species Specificity:

mouse

Background Information

CD28 (T-cell-specific surface glycoprotein CD28), also known as T44 and Tp44, is a 44 kD disulfide-linked homodimeric type I glycoprotein (PMID: 2162180). It is a member of the immunoglobulin superfamily and is expressed on most T lineage cells, NK cell subsets, and plasma cells (PMID: 2162180, 8386518). CD28 may affect in vivo immune responses by functioning both as a cell adhesion molecule linking B and T lymphocytes and as the surface component of a novel signal transduction pathway (PMID: 2162180, 3021470). CD28 binds both CD80 and CD86 with a highly conserved motif MYPPY in the CDR3-like loop (PMID: 15696168, 7964482). CD28 is considered a major co-stimulatory molecule, inducing T lymphocyte activation and IL-2 synthesis, and preventing cell death (PMID: 1348520).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD28 Rabbit Recombinant Antibody (98195-1-RR, Clone: 241707H6) or Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were co-stained with APC Anti-Mouse CD3. Cells were not fixed. This data was developed using the same antibody clone with 98195-1-PBS in a different storage buffer formulation.

1x10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD28 Rabbit Recombinant Antibody (98195-1-RR, Clone:241707H6) (red) or Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue) and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98195-1-PBS in a different storage buffer formulation.