

For Research Use Only

Anti-Human CD207 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98368-2-PBS



Basic Information

Catalog Number:

98368-2-PBS

Concentration:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2333

GenBank Accession Number:

BC022278

GeneID (NCBI):

50489

UNIPROT ID:

Q9UJ71

Full Name:

CD207 molecule, langerin

Calculated MW:

328 aa, 37 kDa

Purification Method:

Protein A purification

CloneNo.:

242303C1

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or GM-CSF, IL-4 and TGF- β treated human PBMCs were intracellularly stained with 0.25 ug Anti-Human CD207 Rabbit RecAb (98368-2-RR, Clone:242303C1) and PE-Conjugated Donkey Anti-Rabbit IgG(H+L). Cells were then stained with APC Anti-Human CD1a (HI 149) (APC-65172, Clone: HI 149). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 98368-2-