

For Research Use Only

Anti-Mouse IL-7Ra/CD127 Rabbit Recombinant Antibody

Catalog Number: 98385-1-RR



Basic Information

Catalog Number:

98385-1-RR

Concentration:

100ug, 1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1826

GenBank Accession Number:

NM_008372.4

GeneID (NCBI):

16197

UNIPROT ID:

P16872

Full Name:

interleukin 7 receptor

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

242852B4

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in 100 µl suspension

Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : mouse splenocytes,

Background Information

CD127, also known as IL-7R subunit alpha (IL-7R α), is a type I membrane glycoprotein expressed on thymocytes, B cell precursors, most T cells, and some lymphoid and myeloid cells (PMID: 8415665). IL-7R is a heterodimer composed of CD127 and IL-2R γ (CD132), which is a common gamma chain shared by the receptors of various cytokines, including interleukins 2, 4, 7, 9, and 15 (PMID: 8266077; 9862091). IL-7R plays critical roles in lymphocyte development and homeostasis (PMID: 26336149). CD127 can also act as a receptor for thymic stromal lymphopoietin (TSLP) (PMID: 10974032).

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse IL-7Ra/CD127 Rabbit RecAb (98385-1-RR, Clone:242852B4) and PE-Conjugated Goat Anti-Rabbit IgG(H+L), or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9). Cells were then stained with APC Anti-Mouse CD3 (17A2) (APC-65077, Clone: 17A2). Cells were not fixed.



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse IL-7Ra/CD127 Rabbit RecAb (98385-1-RR, Clone:242852B4) (red) and PE-Conjugated Goat Anti-Rabbit IgG(H+L), or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue). Cells were not fixed.