

For Research Use Only

Anti-Human CD302 Rabbit Recombinant Antibody

Catalog Number: 98413-1-RR



Basic Information

Catalog Number:

98413-1-RR

Concentration:

100ug, 1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2671

GenBank Accession Number:

NM_014880.5

GeneID (NCBI):

9936

UNIPROT ID:

Q8IX05-1

Full Name:

CD302 molecule

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

242903G5

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.25 ug Anti-Human CD302 Rabbit RecAb (98413-1-RR, Clone: 242903G5) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were then stained with APC Anti-Human CD14. Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



1x10⁶ human PBMCs were surface stained with 0.25 ug Anti-Human CD302 Rabbit RecAb (98413-1-RR, Clone: 242903G5) (red) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.