

For Research Use Only

# Anti-Mouse CD74 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98425-3-PBS



## Basic Information

Catalog Number:

98425-3-PBS

Concentration:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2716

GenBank Accession Number:

NM\_001042605.1

GeneID (NCBI):

16149

UNIPROT ID:

P04441-1

Full Name:

CD74 antigen (invariant polypeptide of major histocompatibility complex, class II antigen-associated)

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

250109G4

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

CD74, also known as MHC class II invariant chain Ii or HLA-DR-associated invariant chain, is a non-polymorphic type II transmembrane glycoprotein. CD74 is mainly distributed in B cells, monocytes, Langerhans cells, dendritic cells and thymic epithelial cells, and is also expressed in epithelial cells. CD74 plays a key role in MHC class II antigen processing and also serves as a cell surface receptor for the macrophage cytokine MIF. CD74 plays an important role in cardiovascular diseases such as atherosclerosis, ischemic heart disease, and myocardial hypertrophy. (PMID: 15475450; PMID: 27752708; PMID: 38992165; PMID: 36712241).

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD74 Rabbit RecAb (98425-3-RR, Clone: 250109G4) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were then stained with CoraLite® Plus 647 Anti-Mouse CD19 (1D3) (CL647-65290, Clone: 1D3). Cells were not fixed. This data was developed using the same antibody clone with 98425-3-PBS in a different storage buffer formulation.



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD74 Rabbit RecAb (98425-3-RR, Clone:250109G4) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98425-3-PBS in a different storage buffer formulation.