

For Research Use Only

Anti-Human LIFR Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98432-2-PBS



Basic Information

Catalog Number: 98432-2-PBS	GenBank Accession Number: NM_001127671.1	Purification Method: Protein A purification
Concentration: 1mg, 2 mg/ml	GeneID (NCBI): 3977	CloneNo.: 250110B3
Source: Rabbit	UNIPROT ID: P42702-1	
Isotype: IgG	Full Name: leukemia inhibitory factor receptor alpha	
Immunogen Catalog Number: EG2599	Calculated MW: 124 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

LIFR, also known as CD118, is a subunit of a receptor for leukemia inhibitory factor (LIF). LIF is a pleiotropic cytokine of the interleukin-6 family that affects the differentiation, survival, and proliferation of a wide variety of cells in the adult and the embryo. LIFR is the low-affinity binding chain that, together with the high-affinity converter subunit gp130, forms a high-affinity receptor complex that mediates the action of LIF. The high-affinity complex also binds a related cytokine, oncostatin M (PMID: 8999038). LIFR has also been identified as a breast cancer metastasis suppressor that functions through the HIPPO-YAP pathway (PMID: 23001183).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 JAR cells were surface stained with 0.25 μ g Anti-Human LIFR Rabbit RecAb (98432-2-RR, Clone: 250110B3) (red) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98432-2-PBS in a different storage buffer formulation.