

For Research Use Only

# Anti-Mouse VISTA Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98487-2-PBS



## Basic Information

Catalog Number:

98487-2-PBS

Concentration:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2809

GenBank Accession Number:

NM\_001159572.1

GeneID (NCBI):

74048

UNIPROT ID:

Q9D659

Full Name:

RIKEN cDNA 4632428N05 gene

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

25012287

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

VISTA, also known as platelet receptor Gi24, B7-H5, C10orf54, Dies1, PD-1H, and SISP1, is a member of the immunoglobulin superfamily. It serves not only as a positive regulator but also as a substrate of MT1-MMP and contributes at least in part to tumor invasion.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug Anti-Mouse VISTA Rabbit RecAb (98487-2-RR, Clone: 250122B7) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9), and PE-conjugated Goat Anti-Rabbit IgG. Cells were then stained with CoraLite® Plus 488 Anti-Mouse CD8a. Cells were not fixed. This data was developed using the same antibody clone with 98487-2-PBS in a different storage buffer formulation.



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug Anti-Mouse VISTA Rabbit RecAb (98487-2-RR, Clone: 250122B7) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-conjugated Goat Anti-Rabbit IgG. Cells were not fixed. This data was developed using the same antibody clone with 98487-2-PBS in a different storage buffer formulation.