

For Research Use Only

Anti-Human TNFRSF21/DR6 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98523-2-PBS



Basic Information

Catalog Number: 98523-2-PBS	GenBank Accession Number: NM_014452.5	Purification Method: Protein A purification
Concentration: 1mg, 3 mg/ml	GeneID (NCBI): 27242	CloneNo.: 250249D8
Source: Rabbit	UNIPROT ID: O75509	
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 21	
Immunogen Catalog Number: EG2718	Calculated MW: 72 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

TNFRSF21 (also known as DR6 or CD358) is an orphan receptor of the tumor necrosis factor (TNF) receptor superfamily. It controls cell death and differentiation in a cell-autonomous manner in different cell types including neurons, oligodendrocytes, endothelial cells, and immune cells (PMID: 29459438). DR6 has been shown to activate NF- κ B and JNK, and induce cell apoptosis (PMID: 9714541). DR6 interacts with TRADD, which is known to mediate signal transduction of TNF receptors. It has been reported that APP binds DR6 to trigger axon pruning and neuron death via distinct caspases, indicating that DR6 might be involved in Alzheimer's disease (PMID: 19225519).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ A431 cells were surface stained with 0.25 ug Anti-Human TNFRSF21/DR6 Rabbit RecAb (98523-2-RR, Clone: 250249D8) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98523-2-PBS in a different storage buffer formulation.



1x10⁶ A431 cells were surface stained with 0.25 ug Anti-Human TNFRSF21/DR6 Rabbit RecAb (98523-2-RR, Clone: 250249D8) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98523-2-PBS in a different storage buffer formulation.