

For Research Use Only

# Anti-Mouse IL-25/IL-17E Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98536-3-PBS



## Basic Information

|  |   |   |
|--|---|---|
| <b>Catalog Number:</b><br>98536-3-PBS      | <b>GenBank Accession Number:</b><br>NM_080729.3 | <b>Purification Method:</b><br>Protein A purification |
| <b>Concentration:</b><br>1mg, 3 mg/ml      | <b>GeneID (NCBI):</b><br>140806                 | <b>CloneNo.:</b><br>250033D3                          |
| <b>Source:</b><br>Rabbit                   | <b>UNIPROT ID:</b><br>Q8VHH8                    |   |
| <b>Isotype:</b><br>IgG                     | <b>Full Name:</b><br>interleukin 25             |   |
| <b>Immunogen Catalog Number:</b><br>EG3316 | <b>Calculated MW:</b><br>19KD                   |   |

## Applications

**Tested Applications:**  
FC (Intra)

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> HEK-293T cells or IL-25/IL-17E transfected HEK-293T cells were intracellularly stained with 0.25 ug Anti-Mouse IL-25/IL-17E Rabbit RecAb (98536-3-RR, Clone: 250033D3) and Multi-rAb CoraLite® Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer. This data was developed using the same antibody clone with 98536-3-PBS in a different storage buffer



1x10<sup>6</sup> HEK-293T cells or IL-25/IL-17E transfected HEK-293T cells were intracellularly stained with 0.25 ug Anti-Mouse IL-25/IL-17E Rabbit RecAb (98536-3-RR, Clone: 250033D3) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and Multi-rAb CoraLite® Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer. This data

