

For Research Use Only

Anti-Human IL12RB2 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98621-2-PBS



Basic Information

Catalog Number: 98621-2-PBS	GenBank Accession Number: NM_001559.2	Purification Method: Protein A purification
Concentration: 1mg, 3 mg/ml	GeneID (NCBI): 3595	CloneNo.: 251356B12
Source: Rabbit	UNIPROT ID: Q99665-1	
Isotype: IgG	Full Name: interleukin 12 receptor, beta 2	
Immunogen Catalog Number: EG4745	Calculated MW: 97 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ PHA-treated human PBMCs were surface stained with 0.25 ug Anti-Human IL12RB2 Rabbit RecAb (98621-2-RR, Clone: 251356B12) (red) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Untreated human PBMCs were surface stained with 0.25 ug Anti-Human IL12RB2 Rabbit RecAb (98621-2-RR, Clone: 251356B12) (black, dashed) and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells



Biolayer interferometry (BLI) kinetic assays of 98621-2-RR against Human IL12RB2 were performed. The affinity constant is 33.7 pM.