

For Research Use Only

Anti-Mouse TGFBI/BIGH3 Rabbit Recombinant Antibody

Catalog Number: 98667-2-RR



Basic Information

Catalog Number:

98667-2-RR

Concentration:

100ug, 1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3160

GenBank Accession Number:

NM_009369.4

GeneID (NCBI):

21810

UNIPROT ID:

P82198

Full Name:

transforming growth factor, beta induced

Calculated MW:

75 kDa

Purification Method:

Protein A purification

CloneNo.:

251191C5

Recommended Dilutions:

FC (Intra): 0.25 ug per 10^6 cells in a 100 μ l suspension

Applications

Tested Applications:

FC (Intra)

Species Specificity:

mouse

Positive Controls:

FC (Intra) : Transfected HEK-293T cells,

Background Information

TGFBI, also named as BIGH3, Kerato-epithelin and RGD-CAP, binds to type I, II, and IV collagens. TGFBI is an adhesion protein which may play an important role in cell-collagen interactions. In cartilage, it may be involved in endochondral bone formation. TGFBI is an extracellular matrix adaptor protein, it has been reported to be differentially expressed in transformed tissues. TGFBI is a predictive factor of the response to chemotherapy, and suggest the use of TGFBI-derived peptides as possible therapeutic adjuvants for the enhancement of responses to chemotherapy. (PMID:20509890) Defects in TGFBI are the cause of epithelial basement membrane corneal dystrophy (EBMD). Defects in TGFBI are the cause of corneal dystrophy Groenouw type 1 (CDGG1). Defects in TGFBI are the cause of corneal dystrophy lattice type 1 (CDL1). Defects in TGFBI are a cause of corneal dystrophy Thiel-Behnke type (CDTB). Defects in TGFBI are the cause of Reis-Buecklers corneal dystrophy (CDRB). Defects in TGFBI are the cause of lattice corneal dystrophy type 3A (CDL3A). Defects in TGFBI are the cause of Avellino corneal dystrophy (ACD).

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ HEK-293T cells (left) or Tgfb1-transfected HEK-293T cells (right) were intracellularly stained with 0.25 ug Anti-Mouse TGFBI/BIGH3 Rabbit RecAb (98667-2-RR, Clone:251191C5), or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) and PE-conjugated Goat Anti-Rabbit IgG. Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).



1x10⁶ HEK-293T cells(left) or Tgfb1-transfected HEK-293T cells(right) were intracellularly stained with 0.25 ug Anti-Mouse Tgfb1 Rabbit RecAb (98667-2-RR, Clone:251191C5) (red), or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue) and PE-conjugated Goat Anti-Rabbit IgG. Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).



Bilayer interferometry (BLI) kinetic assays of 98667-2-RR against Mouse TGFBI/BIGH3 were performed. The affinity constant is 1.52 nM.