

For Research Use Only

Anti-Human Claudin 19 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98681-1-PBS



Basic Information

Catalog Number:

98681-1-PBS

Concentration:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC030524

GeneID (NCBI):

149461

UNIPROT ID:

Q8N6F1

Full Name:

claudin 19

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

251126D7

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 U2OS cells were surface stained with 0.25 ug Anti-Human Claudin 19 Rabbit RecAb (98681-1-RR, Clone: 251126D7) (red) or 0.25 ug Isotype Control (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98681-1-PBS in a different storage buffer formulation.