

For Research Use Only

# Anti-Human Claudin 19 Rabbit Recombinant Antibody

Catalog Number: 98681-1-RR



## Basic Information

Catalog Number:

98681-1-RR

Concentration:

100ug, 1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC030524

GeneID (NCBI):

149461

UNIPROT ID:

Q8N6F1

Full Name:

claudin 19

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

251126D7

Recommended Dilutions:

FC: 0.25 ug per 10<sup>6</sup> cells in a 100 µl suspension

## Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : U2OS cells,

## Background Information

### Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> U2OS cells were surface stained with 0.25 ug Anti-Human Claudin 19 Rabbit RecAb (98681-1-RR, Clone: 251126D7) (red) or 0.25 ug Isotype Control (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed.