

For Research Use Only

# Anti-Human CD179b Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98738-2-PBS



## Basic Information

Catalog Number:

98738-2-PBS

Concentration:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3515

GenBank Accession Number:

BC012293

GeneID (NCBI):

3543

UNIPROT ID:

P15814

Full Name:

immunoglobulin lambda-like  
polypeptide 1

Calculated MW:

213 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

252001D6

## Applications

Tested Applications:

FC

Species Specificity:

human

## Background Information

CD179b, also known as IGLL1, VpreB (surrogate light chain, VpreB1), is a crucial protein component of the pre-B cell receptor (pre-BCR), which is essential for the early development of B lymphocytes (B cells) in the bone marrow.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> NALM-6 cells (left) or REH cells (right) were surface stained with 0.25 ug Anti-Human CD179b Rabbit RecAb (98738-2-RR, Clone: 252001D6) (red) or Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue), and PE-conjugated Goat Anti-Rabbit IgG. Cells were not fixed. This data was developed using the same antibody clone with 98738-2-PBS in a different storage buffer formulation.