

For Research Use Only

Mouse IgG1, λ isotype control Recombinant Antibody

Catalog Number: 98741-1-MR



Basic Information

Catalog Number:

98741-1-MR

Concentration:

100ug, 1000 μ g/ml

Source:

Mouse

Isotype:

IgG1, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

253505A1

Recommended Dilutions:

FC: 0.25 ug per 10^6 cells in a 100 μ l suspension

Applications

Tested Applications:

FC

Species Specificity:

n/a

Positive Controls:

FC : human PBMCs, rat splenocytes

Background Information

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.25 ug Mouse IgG1, λ isotype control RecAb (98741-1-MR, Clone:253505A1) (red in histogram) and PE-conjugated Goat Anti-Mouse IgG, or 0.25 ug Anti-Human CD4 (RPA-T4) (65143-1-Ig, Clone: RPA-T4) (green in histogram), or unstained cells (blue in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.

1x10⁶ rat splenocytes were surface stained with 0.25 ug Mouse IgG1, λ isotype control RecAb (98741-1-MR, Clone:253505A1) (red in histogram) and PE-conjugated Goat Anti-Mouse IgG, or 0.25 ug Anti-Rat CD45 (OX1) (65228-1-Ig, Clone: OX1) (green in histogram), or unstained cells (blue in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.