

For Research Use Only

Anti-Human MUC1/CA15-3 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98901-1-PBS



Basic Information

Catalog Number:

98901-1-PBS

Concentration:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

4582

UNIPROT ID:

P15941

Full Name:

mucin 1, cell surface associated

Purification Method:

Protein A purification

CloneNo.:

261143A1

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or PHA-treated human PBMCs were surface stained with 0.25 ug Anti-Human MUC1/CA15-3 Rabbit RecAb (98901-1-RR, Clone: 261143A1) and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were then stained with APC Anti-Human CD3 (OKT3) Mouse IgG2a Recombinant Antibody (APC-65569, Clone: OKT3). Cells were not fixed. This data was developed using the same antibody clone with 98901-1-PBS in a different storage buffer formulation.



1x10⁶ PHA-treated human PBMCs were surface stained with 0.25 ug Anti-Human MUC1/CA15-3 Rabbit RecAb (98901-1-RR, Clone: 261143A1) (red) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). 1x10⁶ untreated human PBMCs were surface stained with 0.25 ug Anti-Human MUC1/CA15-3 Rabbit RecAb (98901-1-RR, Clone:261143A1) (black, dashed). Cells were not fixed. This data was

