

For Research Use Only

APC Anti-Mouse LPAM-1 (DATK32)

Catalog Number: APC-65101

1 Publications



Basic Information

Catalog Number:

APC-65101

Concentration:

100ug, 200 ug/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC011184

GeneID (NCBI):

16421

UNIPROT ID:

P26011

Full Name:

integrin beta 7

Purification Method:

Affinity purification

CloneNo.:

DATK32

Excitation/Emission maxima
wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

mouse

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qi You	35111148	Front Immunol	FC

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse bone marrow cells were surface stained with 0.25 ug APC Anti-Mouse LPAM-1 (DATK32) (APC-65101, Clone:DATK32) or 0.25 ug APC Rat IgG2a Isotype Control (2A3) (APC-65209, Clone: 2A3), and 0.25 ug Coralite® Plus 488 Anti-Mouse CD45R (B220) (RA3-6B2) (CL488-65139, Clone: RA3-6B2). Cells were not fixed.



1x10⁶ mouse bone marrow cells were surface stained with 0.25 ug APC Anti-Mouse LPAM-1 (DATK32) (APC-65101, Clone: DATK32) or 0.25 ug APC Rat IgG2a Isotype Control (2A3) (APC-65209, Clone: 2A3). Cells were not fixed.