

For Research Use Only

# APC Anti-Mouse LPAM-1 (DATK32)

Catalog Number: APC-65101

1 Publications



## Basic Information

Catalog Number:

APC-65101

Size:

100ug, 0.2 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC011184

GeneID (NCBI):

16421

UNIPROT ID:

P26011

Full Name:

integrin beta 7

Purification Method:

Affinity purification

CloneNo.:

DATK32

Excitation/Emission maxima  
wavelengths:

650 nm / 660 nm

## Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Mouse

Cited Species:

mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Qi You	35111148	Front Immunol	FC

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> C57BL/6 mouse bone marrow cells were surface stained with 0.25 ug FITC-Anti-Human/Mouse CD45R (B220) (FITC-65139, clone RA3-6B2) and 0.25 ug APC-rat IgG2a isotype control. Cells were not fixed.



1X10<sup>6</sup> C57BL/6 mouse bone marrow cells were surface stained with 0.25 ug FITC-Anti-Human/Mouse CD45R (B220) (FITC-65139, clone RA3-6B2) and 0.25 ug APC-Anti-Mouse LPAM-1 (APC-65101, clone DATK32). Cells were not fixed.