

For Research Use Only

APC Anti-Mouse CD24 (M1/69)

Catalog Number: **APC-65121**



Basic Information

Catalog Number:

APC-65121

Size:

100ug, 0.2 mg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

UNIPROT ID:

P24807

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

Excitation/Emission maxima
wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and 0.06 ug APC-Rat IgG2b isotype control (left) or 0.06 ug APC Anti-Mouse CD24 (APC-65121, Clone: M1/69) (right). Cells were not fixed.